

CellTumbler® a low-cost alternative to rock wave bags



Any commercial available wave bag up to 10/20 litre may be used and operated independently with temperature and gas control.




Mix different bags on the same instrument in a design-yourself-kit reducing the instrument complexity. Kits are delivered non-assembled for you to assemble (similar to when you buy furniture from a large Swedish company).

CellCell®



Starter kit




Expansion platform

CellTumbler essential parts are:

- The aluminium T-frame on to which all the parts are bolted
- The platform(s) to which any wave bag may be applied
- The Drive-Unit which handles one or more platforms
- The Heat-Control-Unit with two independent channels for two platforms



Gas control units



Do-it-yourself kit keeps the cost low



Check out the product and pricing at:

www.celltumbler.com

www.cercell.com